

Sigma

Один из мировых лидеров по производству лабораторных центрифуг. Мы предлагаем полный ассортимент лабораторных центрифуг производства компании Sigma: настольные центрифуги, рефрижераторные центрифуги, центрифуги с нагревом, напольные центрифуги, лабораторные центрифуги, медицинские центрифуги, центрифуги малых объемов, центрифуги большой вместимости.

Артикул	Наименование	Цена, Евро, за шт., с НДС
Sigma 1-14, 1-14K		
10014	Sigma 1-14, laboratory table top microcentrifuge, 220-240 V, 50/60 Hz	1322,92
10020	Sigma 1-14K, refrigerated microcentrifuge, 220-240 V, 50/60 Hz	3929,40
Sigma 1-16, 1-16K		
10025	Sigma 1-16, laboratory table top microcentrifuge, 220-240 V, 50/60 Hz	1708,83
1002590	Sigma 1-16, laboratory table top microcentrifuge, incl. rotor 12134, 220-240 V, 50/60 Hz	2122,49
10030	Sigma 1-16K, refrigerated microcentrifuge, 220-240 V, 50/60 Hz	4880,95
1003090	Sigma 1-16K, refrigerated microcentrifuge, incl. rotor 12024, 220-240 V, 50/60 Hz	5781,53
10035	Sigma 1-16, laboratory table top microcentrifuge, 220-240 V, 50/60 Hz, IVD	1802,66
10040	Sigma 1-16K, refrigerated microcentrifuge, 220-240 V, 50/60 Hz, IVD	4953,17
Sigma 2-6, 2-6E		
10208	Sigma 2-6E, laboratory table top centrifuge, 3.900 rpm 220-240 V, 50/60 Hz	2107,95
10223	Sigma 2-6, laboratory table top centrifuge, 4.000 rpm, incl. stainless steel bowl and imbalance switch, 220-240 V, 50/60 Hz	2392,10

Sigma 2-16P, 2-16KL		
10155	Sigma 2-16P , laboratory table top centrifuge, 220-240 V, 50/60 Hz	3752,40
10350	Sigma 2-16KL , refrigerated table top centrifuge, 220-240 V, 50/60 Hz	7162,13
10353	Sigma 2-16KHL , refrigerated table top centrifuge, incl. heating device, max. rotor temp. +40°C up to +60°C depending on rotor and speed, 220-240 V, 50/60 Hz	8503,08

Sigma 3-16L, 3-16KL		
10355	Sigma 3-16L , laboratory table top centrifuge, 220-240 V, 50/60 Hz	5363,81
10360	Sigma 3-16KL , refrigerated table top centrifuge, 220-240 V, 50/60 Hz	7398,60
10358	Sigma 3-16L , laboratory table top centrifuge, 220-240 V, 50/60 Hz, IVD	5464,34
10363	Sigma 3-16KL , refrigerated table top centrifuge, 220-240 V, 50/60 Hz, IVD	7499,14

Sigma 3-18KS		
10370	Sigma 3-18KS , refrigerated table top centrifuge, 220-240 V, 50/60 Hz	9089,30
10374	Sigma 3-18KHS , refrigerated table top centrifuge, incl. heating device, max. rotor temperature +40°C up to +60°C depending on rotor and speed, 220-240 V, 50/60 Hz	10431,67

Sigma 3-30KS		
10375	Sigma 3-30KS , refrigerated high speed table top centrifuge, 220-240 V, 50 Hz	13004,54
10379	Sigma 3-30KHS , refrigerated high speed table top centrifuge, incl. heating device, max. rotor temp. +40°C up to +60°C depending on rotor and speed, 220-240 V, 50 Hz	14436,40

Sigma 4-5L, 4-16S, 4-16KS		
10405	Sigma 4-5L , laboratory table top centrifuge, 220-240 V, 50/60 Hz	5857,99
1040590	Sigma 4-5L , laboratory table top centrifuge, incl. rotor no. 11650, 2 sets buckets no. 13450, 2 sets caps no. 17170, 2 sets adapter no. 17658, 220-240 V, 50/60 Hz	9893,59
10407	Sigma 4-5L , laboratory table top centrifuge, 220-240 V, 50/60 Hz, IVD	5930,21
10380	Sigma 4-16S , laboratory table top centrifuge, 220-240 V, 50/60 Hz	7094,16
10385	Sigma 4-16KS , refrigerated table top centrifuge, 220-240 V, 50 Hz	10900,56
10389	Sigma 4-16KHS , refrigerated table top centrifuge, incl. heating device, max. rotor temperature +40°C up to +60°C depending on rotor and speed, 220-240 V, 50 Hz	12770,43

Sigma 6-16S, 6-16KS		
10390	Sigma 6-16S , laboratory table top centrifuge, 220-240 V, 50/60 Hz	9794,94
10393	Sigma 6-16HS , laboratory table top centrifuge, incl. heating device up to +60°C; max.+70°C with rotor no. 11165, 220-240 V, 50/60 Hz	11972,94
10395	Sigma 6-16KS , refrigerated table top centrifuge, 220-240 V, 50 Hz	14884,99
10399	Sigma 6-16KHS , refrigerated table top centrifuge, incl. heating device, max. rotor temp. +40°C up to +60°C depending on rotor and speed, 220-240 V, 50 Hz	16495,08

Sigma 8KS		
10621	Sigma 8KS , floorstanding refrigerated centrifuge, aircooled compressor, 3 x 400 V, 50 Hz	40094,32
10623	Sigma 8KS GMP , floorstanding refrigerated centrifuge, watercooled compressor and electronic, gasketed casing, 3 x 400 V, 50 Hz	43032,71
10944	Sigma 8KBS no. 10625, floorstanding refrigerated centrifuge, aircooled compressor, 3 x 400 V, 50 Hz , incl. rotor no. 11805, 3 sets buckets no. 13860 with 3 sets adapters no. 13865, interface RS 232	48634,88

Роторы затухания		
-------------------------	--	--

11024	Microhematocritrotor incl. lid, suitable for 24 capillary tubes Ø 1.5 x 75 mm, 50 µl no. 15001, incl. reader no. 17029, usable with reader no. 17024	658,44
11026	Microhematocritrotor incl. lid, suitable for 24 capillary tubes Ø 1.3 x 50 mm, 19 µl no. 15028, incl. reader no. 17028, usable with reader no. 17026	291,70
11030	Swing-out rotor, 4 place, for buckets no. 13031, 13033, 13034, 13035, 13037, 13038, 13099; buckets no. 13034, 13037 for rotors from the year 2003 onwards only	553,66
11031	Swing-out rotor 16 x 15 ml complete, no. 11030, incl. 2 sets carriers no. 13035 and 16 PS-tubes 15 ml no. 15020	807,12
11032	Swing-out rotor 4 x 100 ml complete, no. 11030, incl. 2 sets buckets no. 13099 and 4 glass tubes 100 ml no. 15100	907,66
11035	Swing-out rotor 8 x 15 ml complete, no. 11036, incl. 1 set carriers no. 13035 and 8 PS- tubes 15 ml no. 15020	288,86
11036	Swing-out rotor, 2 place, for carriers no. 13031, 13033, 13034, 13035, 13038, 13099	189,42
11070	Swing-out rotor, 4 place, for bucket no. 13299, Software version 22 or higher	553,66
11118	Swing-out rotor for microtiter plates, Biosafe certificate, incl. 1 set carriers no. 13218, with 1 set caps no. 17108, incl. plate holder no. 17978, max. height of plates 80 mm, max. load per bucket approx. 520 grams	3098,21

11121	Swing-out rotor for microtiter plates, incl. 1 set carriers no. 13223, max. height of plates 50 mm, max. load per bucket approx. 345 grams	1818,14
11122	Swing-out rotor for microtiter plates, incl. 1 set carriers no. 13222, max. height of plates 56 mm, max. load per bucket approx. 335 gram	1533,53
11123	Swing-out rotor for microtiter plates, incl. 1 set carriers no. 13223, max. height of plates 50 mm, max. load per bucket approx. 345 gram	1907,35
11124	Swing-out rotor 24 x 1.5 / 2.0 ml, incl. 3 sets buckets no. 13124, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17878	1635,48
11128	Swing-out rotor 6 x 1.5 / 2.0 ml, incl. 3 sets buckets no. 13128, for e.g. reaction vials no. 15008, 15040	1144,13
11132	Swing-out rotor 28 x 15 ml complete, no. 11133, incl. 2 sets buckets no. 13104, 2 sets round carriers no. 17015 and 28 PP-tubes 15 ml no. 15023	1416,00
11133	Swing-out rotor, 4 place, for buckets no. 13104, 13130	659,86
11134	Swing-out rotor, 4 x 10 ml, incl. 2 sets buckets no. 13134 for e.g. tubes no. 15000, 15039	2098,51
11135	Swing-out rotor 24 x 1.5/2.0 ml, incl. 3 sets buckets no. 13124, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17879	1743,10
11136	Drum rotor for 6 cassettes no. 14000, 14002, 14020, 14021, for reaction vials (max. 60/192 pcs.), incl. hermetic aluminium lid no. 17898	1231,92
11148	Swing-out rotor 32 x 1.5/2.0 ml, incl. 4 sets buckets no. 13124, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17880	2439,77
11150	Swing-out rotor, 4 place, for buckets no. 13215, 13220, 13221, 13233, 13234, 13235, 13236, 13350, (13420, 13550 with 6-16S, 6-16KS only)	1141,30
11151	Swing-out rotor 48 x 15 ml complete, no. 11150, incl. 2 sets rectangular buckets no. 13215, 2 sets rectangular carriers no. 18015 and 48 PPCO-tubes 15 ml no. 15023	2334,98
11152	Swing-out rotor 96 x 15 ml complete, no. 11150, incl. 1 set double rectangular buckets no. 13420, 4 sets rectangular carriers no. 18015 and 96 PPCO-tubes 15 ml no. 15023	3790,63
11156	Swing-out rotor, 6 place, for buckets no. 13127	1506,62
11157	Swing-out rotor 72 x 15 ml complete, no. 11156, incl. 3 sets rectangular buckets no. 13127, 3 sets rectangular carriers no. 18015 and 72 PPCO-tubes 15 ml no. 15023	2843,33
11165	Swing-out rotor, 4 place, for buckets no. 13113, 13114, 13116 for glass tubes for mineral oil	1175,28
11166	Swing-out rotor 4 x 100 ml complete, no. 11165, incl. 2 sets buckets no. 13113, 4 mineral oil glass tubes 8" following former ASTM D91/96 no. 15130 and 2 sets adapters no. 16103, (13113 sealable with cap no. 17107)	2806,51
11170	Swing-out rotor, 4 place, for bucket no. 13299, 2-16P software 27 or higher; 2-16KL software 32 or higher	801,46

11171	Swing-out rotor, 4 place, for bucket no. 13299, 3-30KS: Software 47 or higher	955,80
11180	Swing-out rotor, 4 place, for buckets no. 13180, 13190, 13194	801,46
11181	Swing-out rotor 48 x 15 ml complete, no. 11180, incl. 2 sets rectangular buckets no. 13180, 2 sets rectangular carriers no. 18015 and 48 PPCO-tubes 15 ml no. 15023	1969,66
11182	Swing-out rotor 12 x 50 ml culture tubes complete, no. 11180, incl. 2 sets buckets no. 13190, 2 sets adapters no. 17346 and 12 tubes no. 15151	2478,00
11190	Swing-out rotor 4 x 100 ml complete, no. 11192, incl. 2 sets buckets no. 13097 and 4 glass tubes 100 ml no. 15100	1360,78
11191	Swing-out rotor 16 x 15 ml complete, no. 11192, incl. 2 sets carriers no. 13012 and 16 PS-tubes 15 ml no. 15020	1397,59
11192	Swing-out rotor, 4 place, for buckets no. 13009, 13012, 13041, 13042, 13097, 13150, 13152; buckets 13041 and 13152 for rotors from the year 2003 onwards only	801,46
11194	Swing-out cytology rotor complete, incl. 2 sets buckets no. 13224, 4 inserts no. 15223	3091,13
11222	Swing-out rotor for microtiter plates, incl. 1 set carriers no. 13222, max. height of plates 56 mm, max. load per bucket approx. 335 grams	1505,21
11224	Swing-out cytology rotor complete, incl. 2 sets buckets no. 13224, 4 inserts no. 15223	3140,69
11234	Swing-out cytology rotor complete, incl. 2 sets buckets no. 13224, 4 inserts no. 15223	3140,69
11241	Swing-out rotor for microtiter plates, Biosafe certificate, incl. 1 set buckets no. 13241, 1 set caps no. 17108 and plate holder no. 17982, max. height of plates 55 mm, without cap max. 66 mm, max. load per bucket approx. 290 grams	1870,54
11390	Swing-out rotor, 4 place, for buckets, no. 13009, 13012, 13041, 13042, 13097, 13150, 13152; buckets 13041 and 13152 for rotors from the year 2003 onwards only	955,80
11391	Swing-out rotor 4 x 100 ml complete, no. 11390, incl. 2 sets buckets no. 13097 and 4 glass tubes 100 ml no. 15100	1515,12
11392	Swing-out rotor 16 x 15 ml complete, no. 11390, incl. 2 sets carriers no. 13012 and 16 PS-tubes 15 ml no. 15020	1540,61
11409	Microhematocritrotor suitable for 24 capillary tubes Ø 1.5 x 75 mm, 50 µl no. 15001, incl. reader no. 17029, usable with reader 17002, incl. rotor screw 81543	1331,04
11650	Swing-out rotor, 4 place, for buckets 13421, 13450, (13621, 13625, 13635, 13650 only in 6-16S, 6-16KS)	1689,29
11651	Swing-out rotor 48 x 15 ml culture tubes complete, no. 11650, incl. 2 sets buckets no. 13650, 2 sets adapters no. 17659 and 48 tubes 15 ml no. 15115	4146,05
11660	Swing-out rotor, stainless steel, 4 place, for buckets 13421, 13450, (13460 only in 4-5L), (13621, 13625, 13635, 13650 only in 6-16S, 6-16KS)	1689,29
11805	Swing-out rotor with windshield, 6 place, for buckets no. 13845, 13850, 13860	11451,85

11806	Swing-out rotor, 6 place, for buckets no. 13845, 13850, 13855, 13860	8425,20
-------	----------------------------------------------------------------------	----------------

Угловые роторы

12011	Angle rotor 6 x 15 ml complete, incl. 3 sets buckets no. 13011, for tubes Ø 17 x up to 130 mm, e.g. no. 15015, 15020, 15023, 15024, culture tubes 15 ml no. 15115 and for 6 reaction vials 1.5/2.0 ml no. 15008, 15040	529,58
12012	Angle rotor 8 x 15 ml complete, incl. 4 sets buckets no. 13011, for tubes Ø 17 x up to 115 mm, e.g. no. 15015, 15020, 15023, 15024 and for 8 reaction vials 1.5/2.0 ml no. 15008, 15040	564,98
12024	Angle rotor 24 x 1.5/2.0 ml, for reaction vials e.g. no. 15008, 15040, incl. polysulfone lid no. 17828	906,24
12061	Angle rotor 30 x 15 ml complete, incl. 15 sets buckets no. 13011, for tubes Ø 17 x up to 120 mm, e.g. no. 15015, 15020, 15023, 15024, Monovettes and culture tubes 15 ml no. 15115	924,65
12062	Angle rotor 20 x 15 ml complete, incl. 10 sets buckets no. 13011, for tubes Ø 17 x up to 120 mm, e.g. no. 15015, 15020, 15023, 15024, Monovettes and culture tubes 15 ml no. 15115	865,18
12071	Angle rotor for 12 x 15 ml culture tubes, e.g. no. 15115 and 6 x 50 ml culture tubes, e.g. no. 15151, incl 6 sets buckets 13011, 3 sets 13151	1525,03
12072	Angle rotor 20 x 15 ml complete, incl. 10 sets buckets no. 13011, for tubes Ø 17 x up to 120 mm, e.g. no. 15015, 15020, 15023, 15024, Monovettes and culture tubes 15 ml no. 15115	918,98
12073	Angle rotor 30 x 15 ml complete, incl. 15 sets buckets no. 13011, for tubes Ø 17 x up to 120 mm, e.g. no. 15015, 15020, 15023, 15024, Monovettes and culture tubes 15 ml no. 15115	975,62
12082	Angle rotor 12 x 1.5/2.0 ml, for reaction vials e.g. no. 15008, 15040, incl. polysulfone lid no. 17882	376,66
12083	Angle rotor 18 x 1.5/2.0 ml, for reaction vials e.g. no. 15008, 15040 or 18 x 0.25/0.4 ml e.g. no. 15014, 2 lines, incl. polysulfone lid no. 17882	412,06
12084	Angle rotor 24 x 1.5/2.0 ml, for reaction vials e.g. no. 15008, 15040, incl. polysulfone lid no. 17882	412,06
12092	Angle rotor made of polypropylene, 12 x 1.5/2.0 ml, for reaction vials e.g. no. 15008, 15040, incl. polysulfone lid no. 17882	211,55
12093	Angle rotor made of polypropylene, 18 x 1.5/2.0 ml, for reaction vials e.g. no. 15008, 15040 or 18 x 0,25/0,4 ml e.g. no. 15014, 2 lines, incl. polysulfone lid no. 17882	266,63
12094	Angle rotor made of polypropylene, 24 x 1.5/2.0 ml, for reaction vials e.g. no. 15008, 15040, incl. polysulfone lid no. 17882	266,63

12096	Angle rotor made of polypropylene, for 2 PCR-strips with 8 tubes 0.2 ml e.g no. 15042, incl. polysulfone lid no. 17882	300,19
12104	Angle rotor, for 12 PCR-strips with 8 tubes 0.2 ml e.g no. 15042, incl. hermetic aluminium lid no. 17867	1230,50
12110	Angle rotor 12 x 1.5/2.0 ml, for e.g. reaction vials e.g. no. 15008, 15040, incl. hermetic aluminium lid no. 17834	790,13
12111	Angle rotor 10 x 10 ml, for tubes no. 15000, 15010, 15039, incl. hermetic aluminium lid no. 17835	1134,22
12115	Angle rotor, for 12 PCR-strips with 8 tubes 0.2 ml e.g no. 15042, incl. hermetic aluminium lid no. 17875	1982,40
12130	Angle rotor 30 x 1.5/2.0 ml, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17844	1315,46
12131	Angle rotor 30 x 1.5/2.0 ml, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17848	1315,46
12132	Angle rotor 30 x 1.5/2.0 ml, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17849	1078,99
12133	Angle rotor 30 x 1.5/2.0 ml, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17863	1078,99
12134	Angle rotor made of polypropylene 24 x 1.5/2.0 ml, for e.g. reaction vials no. 15008, 15040, incl. polysulfone lid no. 17929	652,78
12135	Angle rotor made of polypropylene, for 4 PCR-strips with 8 tubes 0.2 ml e.g no. 15042, incl. polysulfone lid no. 17930	948,72
12136	Angle rotor 36 x 1.5/2.0 ml, for reaction vials e.g. no. 15008, 15040, incl. lid no. 17613	1323,96
12137	Angle rotor made of polypropylene 24 x 1.5/2.0 ml, for e.g. reaction vials no. 15008, 15040, incl. polysulfone lid no. 17932, extra cool	948,72
12138	Angle rotor 4 x 80/85 ml, for tubes e.g. no. 15076, 15080, incl. hermetic aluminium lid no. 17897	1335,29
12139	Angle rotor 6 x 30 ml, for tubes no. 15029, 15030, 15032, incl. hermetic aluminium lid no. 17851	1151,21
12141	Angle rotor 10 x 10 ml, for tubes no. 15000, 15010, 15039, incl. hermetic aluminium lid no. 17850	1043,59
12148	Angle rotor 24 x 1.5/2.0 ml, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17864	938,81
12150	Angle rotor 6 x 50 ml, for tubes e.g. no. 13055, 15051, 15052, incl. hermetic aluminium lid no. 17837	1525,03
12154	Angle rotor 24 x 1.5/2.0 ml, for e.g. reaction vials no. 15008, 15040, incl. hermetic aluminium lid no. 17829	938,81
90934	Angle rotor 24 x 1.5/2.0 ml, similar rotor 12154, Cleanroom coating, incl. coated hermetic aluminium lid no. 90933	1220,59
12155	Angle rotor 4 x 80/85 ml, for tubes e.g. no. 13085, 15076, 15080, incl. hermetic aluminium lid no. 17839	1533,53
12156	Angle rotor 8 x 50 ml, for tubes e.g. no. 13055, 15051, 15052, incl. hermetic aluminium lid no. 17840	1533,53
12157	Angle rotor 20 x 10 ml, for tubes e.g. no. 15000, 15010, 15039, incl. hermetic aluminium lid no. 17841	1533,53
12158	Angle rotor 6 x 30 ml, for tubes e.g. no. 15029, 15030, 15032, incl. hermetic aluminium lid no. 17842	1151,21
12159	Angle rotor 6 x 85 ml, for tubes e.g. no. 13085, 15076, 15080, incl. hermetic aluminium lid no. 17843	1764,34

12168	Angle rotor 20 x 10 ml, for tubes no. 15000, 15010, 15039, incl. hermetic aluminium lid no. 17846	2033,38
12170	Angle rotor for 12 x 15 ml culture tubes e.g. no. 15115, incl. hermetic aluminium lid no. 17859	2143,82
12171	Angle rotor for 12 x 15 ml culture tubes e.g. no. 15115, incl. hermetic aluminium lid no. 17872	1561,85
12172	Angle rotor 12 x 30 ml, for tubes no. 15029, 15030, 15032, incl. hermetic aluminium lid no. 17845	1852,13
12174	Angle rotor for 24 x 15 ml culture tubes e.g. no. 15115, incl. hermetic aluminium lid no. 17883	2410,03
12180	Angle rotor 12 x 5 ml, Eppendorf reaction caps, incl. hermetic aluminium lid no. 17180	938,81
12181	Angle rotor for 6 x 50 ml culture tubes, e.g no. 15151, incl. hermetic aluminium lid no. 17931	1319,71
12182	Angle rotor 12 x 5 ml, Eppendorf reaction caps, incl. hermetic aluminium lid no. 17182 (3-16L, 3-16KL Software version 32 or higher, 3-18KS Software version 46 or higher)	938,81
12201	Angle rotor 100 x 1.5/2.0 ml, for e.g. reaction vials e.g. no. 15008, 15040 (2 lines, angle 45°), incl. hermetic aluminium lid no. 17896	3300,70
12265	Angle rotor 6 x 85 ml, Biosafe with certificate, for tubes no. 13085, 15076, 15080, incl. bayonet lid no. 17607	1728,94
12266	Angle rotor 8 x 50 ml, Biosafe with certificate, for tubes no. 13055, 15051, 15052, incl. bayonet lid no. 17605	1554,77
12269	Angle rotor for 8 x 50 ml, Biosafe with certificate, culture tubes e.g. no. 15151, incl. bayonet lid no. 17604	2033,38
12349	Angle rotor 48 x 1.5/2.0 ml, for reaction vials, e.g. no. 15008, 15040, incl. hermetic lid 17349	1557,60
12356	Angle rotor 6 x 250 ml, Biosafe with certificate, for bottles no. 13255, 15247, 15249 and adapters 14122 - 14134, incl. bayonet lid no. 17603	3649,03
12449	Angle rotor 48 x 1.5/2.0 ml, Biosafe with certificate, for reaction vials no. e.g. 15008, 15040, incl. bayonet lid no. 17612	1687,87
12450	Angle rotor 24 x 5 ml Eppendorf reaction caps, Biosafe with certificate, incl. bayonet lid no. 17614	1687,87
12505	Angle rotor 6 x 500 ml, for bottles no. 13507, 15508 and adapters 14144 - 14155, sealable with hermetic aluminium lid no. 17891	6072,56
12510	Angle rotor 6 x 1000 ml, incl. steel buckets no. 13511, for bottles no. 13510, 15921 and adapters 14162 - 14182, sealable with hermetic aluminium lid no. 17892	11058,68
12600	Angle rotor 6 x 500 ml, Biosafe with certificate, for bottles no. 13507, 15508 and adapters 14144 - 14155, incl. bayonet lid no. 17600, incl. torque wrench set 17060	6200,66
90938	Angle rotor 6 x 500 ml, similar rotor 12600, Cleanroom coating, incl. coated bayonet lid no. 90937, incl. torque wrench set 17060	8059,87

19776	Angle rotor for 6 x 50 ml culture tubes e.g. no. 15151, incl. hermetic aluminium lid no. 17831	1525,03
90940	Angle rotor 6 x 50 ml culture tubes, similar rotor 19776, Cleanroom coating, incl. coated hermetic aluminium lid no. 90939	1982,40